

Hypoxic Culture of HeLa Cells Using MakCell

The tabletop automated culture system MakCell is an automated medium exchange device that also supports hypoxic culture. In conventional hypoxic culture, when medium is exchanged manually, cultured cells must be exposed to ambient air even if the incubator itself is under hypoxic conditions. In contrast, MakCell performs medium exchange while maintaining the hypoxic environment, enabling stable cell culture without disturbing the culture conditions (Ref. Application Note "[Stability in hypoxic environments of MakCell](#)"). In this study, HeLa cell lines were cultured and analyzed using the hypoxia function of MakCell. Immunofluorescence staining was performed with an antibody against HIF1 α , a protein known to be specifically localized in the nucleus under hypoxic conditions. As a result, nuclear-specific localization of HIF1 α was observed starting at approximately 1.0 hour of culture.

METHOD

- Cell type: HeLa cell line (seeded at 2×10^5 cells/well after maintenance culture in T-75 flasks, followed by trypsinization and seeding into the culture vessels)
- Maintenance medium: DMEM (high glucose, Wako), 10% FBS (Sigma), 1% Anti-Anti (Gibco)
- Culture vessels: 12-well plate (Falcon), ϕ 15 glass (Matsunami)
- Cell counting device: ViCELL MetaFLEX (Beckman Coulter)
- Antibodies: HIF1 α antibody (Abcam), fluorescently labeled secondary antibody (ThermoFisher)


Cells were transferred into MakCell set to 1% O₂, then fixed and fluorescently stained at 1, 3, and 6 h.

RESULTS


- Measurement of medium inside the incubator conditions:

MakCell culture: CO₂ 5%, O₂ 1%

Standard culture: CO₂ 5%, O₂ 19–20%



- Localization of HIF1 α in HeLa cells under each oxygen concentration




HIF1 α
DAPI

<Reference> Localization of HIF1 α in human adipose-derived mesenchymal stromal cells under each oxygen concentration

Conditions:

- MakCell culture: CO₂ 5%, O₂ 1%
- Standard culture: CO₂ 5%, O₂ 19–20%

Samples: 3-hour samples under both conditions, using the same culture vessels and equipment as above.



DISCUSSION

- MakCell was set to a hypoxic environment, achieving a hypoxic condition in the medium within just 1 hour and subsequently maintaining a stable state.
- Immunofluorescence staining confirmed that HIF1 α translocates to and localizes within the nucleus under a hypoxic condition.
- MakCell demonstrated that it can accommodate experiments measuring the temporal response of HIF1 α nuclear localization during transition to a hypoxic environment. This is because, after reaching the set O₂ concentration, fluctuations in O₂ concentration caused by opening and closing the door during storage of the cell culture vessel inside the incubator return to the set value within a relatively short time of 12 minutes (Ref. Application Note "[Stability in hypoxic environments of MakCell](#)").
- Long-term culture under hypoxic conditions using MakCell offers the advantage of avoiding O₂ concentration fluctuations and maintaining stability in the cell culture environment, as the automatic medium exchange function eliminates the need to take the culture vessel out from the incubator. We plan to evaluate the usefulness of MakCell for automated medium exchange under hypoxic conditions, including the use of mesenchymal stromal cells, which have been reported to exhibit enhanced proliferation and stem cell marker gene expression under hypoxia.

